

ALAS1 Monoclonal Antibody



Catalog Number: 160696

Product name

ALAS1 Monoclonal Antibody

Specificity

Human

Antibody description

Monoclonal antibody to ALAS1

Preparation

Antigen: Recombinant protein of human ALAS1

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

Validations




Western blot - ALAS1 Monoclonal Antibody

Western blot analysis of extracts of various cell lines, using ALAS1 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - ALAS1 Monoclonal Antibody

Immunofluorescence analysis of MCF-7 cells using ALAS1 antibody.



Immunofluorescence - ALAS1 Monoclonal


ALAS1 Monoclonal Antibody



Catalog Number: 160696

Antibody

Immunofluorescence analysis of HeLa cells using ALAS1 antibody .



Immunofluorescence - ALAS1 Monoclonal Antibody

Immunofluorescence analysis of HepG2 cells using ALAS1 antibody .