

KIF23 Monoclonal Antibody



Catalog Number: 160719

Product name

KIF23 Monoclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Monoclonal antibody to KIF23

Preparation

Antigen: Recombinant protein of human KIF23

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

Validations




Western blot - KIF23 Monoclonal Antibody

Western blot analysis of extracts of various cell lines, using KIF23 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - KIF23 Monoclonal Antibody

Immunofluorescence analysis of HeLa cells using KIF23 antibody .



Immunofluorescence - KIF23 Monoclonal Antibody

Immunofluorescence analysis of A549 cells using

KIF23 Monoclonal Antibody



Catalog Number: 160719

KIF23 antibody .