

TG Polyclonal Antibody



Catalog Number: 160813

Product name

TG Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to TG

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 2469-2768 of human TG (NP_003226.4).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - TG Polyclonal Antibody

Western blot analysis of extracts of HeLa cells, using TG antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.