

JUND Polyclonal Antibody



Catalog Number: 160897

Product name

Rabbit IgG

JUND Polyclonal Antibody

Applications

Specificity

WB, IHC

Human, Mouse, Rat

Dilutions

Antibody description

WB 1:500 - 1:2000


IHC 1:50 - 1:100

Polyclonal antibody to JUND

Validations

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human JUND (NP_005345.3).



Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Western blot - JUND Polyclonal Antibody

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Western blot analysis of extracts of various cell lines, using JUND antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.

Clonality

Polyclonal

Ig Type