

CCL26 Polyclonal Antibody



Catalog Number: 160906

Product name

CCL26 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to CCL26

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 24-94 of human CCL26 (NP_006063.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - CCL26 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using CCL26 antibody . Blue: DAPI for nuclear staining.