

TIMP3 Polyclonal Antibody



Catalog Number: 160910

Product name

TIMP3 Polyclonal Antibody

Specificity

Human, Rat

Antibody description

Polyclonal antibody to TIMP3

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 100 to the C-terminus of human TIMP3 (NP_000353.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, FC

Dilutions

WB 1:500 - 1:2000

FC 1:50 - 1:200

Validations



Western blot - TIMP3 Polyclonal Antibody

Western blot analysis of extracts of 293T cells, using TIMP3 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.