

GSR Polyclonal Antibody



Catalog Number: 160944

Product name

GSR Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to GSR

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 273-522 of human GSR (NP_000628.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:200

Validations



Western blot - GSR Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using GSR antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.