

# DEFB121 Polyclonal Antibody



Catalog Number: 160945

## Product name

DEFB121 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to DEFB121

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-76 of human DEFB121 (NP\_001011878.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


IHC, IF

## Dilutions

IHC 1:50 - 1:100


IF 1:50 - 1:100

## Validations




Immunohistochemistry - DEFB121 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human prostate using DEFB121 antibody at dilution of 1:100 (40x lens).



Immunofluorescence - DEFB121 Polyclonal Antibody

Immunofluorescence analysis of rat testis cells using DEFB121 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence - DEFB121 Polyclonal

# DEFB121 Polyclonal Antibody



Catalog Number: 160945

---

Antibody

using DEFB121 antibody at dilution of 1:100.  
Blue: DAPI for nuclear staining.

Immunofluorescence analysis of mouse testis cells