

DEFB121 Polyclonal Antibody



Catalog Number: 160945

Product name

DEFB121 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to DEFB121

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-76 of human DEFB121 (NP_001011878.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC, IF

Dilutions

IHC 1:50 - 1:100


IF 1:50 - 1:100

Validations




Immunohistochemistry - DEFB121 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human prostate using DEFB121 antibody at dilution of 1:100 (40x lens).



Immunofluorescence - DEFB121 Polyclonal Antibody

Immunofluorescence analysis of rat testis cells using DEFB121 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence - DEFB121 Polyclonal

DEFB121 Polyclonal Antibody



Catalog Number: 160945

Antibody

using DEFB121 antibody at dilution of 1:100.

Blue: DAPI for nuclear staining.

Immunofluorescence analysis of mouse testis cells