

TSPYL1 Polyclonal Antibody



Catalog Number: 160981

Product name

TSPYL1 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to TSPYL1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human TSPYL1 (NP_003300.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB 1:500 - 1:2000


IP 1:50 - 1:100

Validations



Western blot - TSPYL1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TSPYL1 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .Exposure time: 3s.

Immunoprecipitation - TSPYL1 Polyclonal Antibody

Immunoprecipitation analysis of 200ug extracts of HeLa cells, using 3 ug TSPYL1 antibody . Western blot was performed from the immunoprecipitate using TSPYL1 antibody at a dilution of 1:1000.