

ACD Polyclonal Antibody



Catalog Number: 160983

Product name

Rabbit IgG

ACD Polyclonal Antibody

Applications

Specificity

WB

Human, Mouse, Rat

Dilutions

Antibody description


WB 1:500 - 1:2000

Polyclonal antibody to ACD

Validations

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human ACD (NP_001075955.1).



Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Western blot - ACD Polyclonal Antibody

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Western blot analysis of extracts of various cell lines, using ACD antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.

Clonality

Polyclonal

Ig Type