

KCNK4 Polyclonal Antibody



Catalog Number: 161004

Product name

KCNK4 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to KCNK4

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 25-90 of human KCNK4 (NP_201567.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - KCNK4 Polyclonal Antibody

Western blot analysis of extracts of HepG2 cells, using KCNK4 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.