

CLCN7 Polyclonal Antibody



Catalog Number: 161071

Product name

CLCN7 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to CLCN7

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 626-805 of human CLCN7 (NP_001278.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:100

Validations



Immunofluorescence - CLCN7 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using CLCN7 antibody . Blue: DAPI for nuclear staining.