

HEXA Polyclonal Antibody



Catalog Number: 161081

Product name

HEXA Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to HEXA

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 270-529 of human HEXA (NP_000511.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Immunohistochemistry - HEXA Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat testis using HEXA antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - HEXA Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using HEXA antibody at dilution of 1:100



(40x lens).


Immunohistochemistry - HEXA Polyclonal Antibody

HEXA Polyclonal Antibody



Catalog Number: 161081

Immunohistochemistry of paraffin-embedded mouse testis using HEXA antibody at dilution of




1:100 (40x lens).

Immunohistochemistry - HEXA Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse spinal cord using HEXA antibody at

dilution of 1:100 (40x lens).



Immunohistochemistry - HEXA Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using HEXA antibody at dilution of 1:100 (40x lens).