

NDP Polyclonal Antibody



Catalog Number: 161093

Product name

NDP Polyclonal Antibody

Specificity

Rat

Antibody description

Polyclonal antibody to NDP

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 25-133 of human NDP (NP_000257.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions

IHC 1:50 - 1:100

Validations



Immunohistochemistry - NDP Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using NDP antibody at dilution of 1:100 (40x lens).