

NTF3 Polyclonal Antibody



Catalog Number: 161094

Product name

NTF3 Polyclonal Antibody

Specificity

Mouse, Rat

Antibody description

Polyclonal antibody to NTF3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 19-257 of human NTF3 (NP_002518.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:100

Validations



Immunohistochemistry - NTF3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat lung using NTF3 antibody at dilution of 1:100 (40x



lens).

Immunohistochemistry - NTF3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using NTF3 antibody at dilution of 1:100 (40x lens).