

# RAD51C Polyclonal Antibody



Catalog Number: 161100

---

## Product name

RAD51C Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to RAD51C

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-135 of human RAD51C (NP\_002867.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions

IF 1:50 - 1:100

## Validations



Immunofluorescence - RAD51C Polyclonal Antibody

Immunofluorescence analysis of MCF7 cells using RAD51C antibody .