

RPN1 Polyclonal Antibody



Catalog Number: 161105

Product name

RPN1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to RPN1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 20-280 of human RPN1 (NP_002941.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Western blot - RPN1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using RPN1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - RPN1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon using RPN1 antibody at dilution of



1:100 (40x lens).


RPN1 Polyclonal Antibody



Catalog Number: 161105


Immunohistochemistry - RPN1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using RPN1 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - RPN1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat kidney using RPN1 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - RPN1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human tonsil using RPN1 antibody at dilution of



1:100 (40x lens).


Immunohistochemistry - RPN1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver cancer using RPN1 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - RPN1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using RPN1 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - RPN1 Polyclonal


RPN1 Polyclonal Antibody



Catalog Number: 161105

Antibody

Immunohistochemistry of paraffin-embedded mouse testis using RPN1 antibody at dilution of



1:100 (20x lens).

Immunohistochemistry - RPN1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse spinal cord using RPN1 antibody at dilution of 1:100 (40x lens).