

SLC30A1 Polyclonal Antibody



Catalog Number: 161117

Product name

SLC30A1 Polyclonal Antibody

Specificity

Mouse, Rat

Antibody description

Polyclonal antibody to SLC30A1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 328-507 of human SLC30A1 (NP_067017.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - SLC30A1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using SLC30A1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.