

# ATF6 Polyclonal Antibody



Catalog Number: 161142

## Product name

ATF6 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to ATF6

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 570 to the C-terminus of human ATF6 (NP\_031374.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

## Validations




Western blot - ATF6 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ATF6 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.



Immunohistochemistry - ATF6 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human esophagus using ATF6 antibody (40x lens).



# ATF6 Polyclonal Antibody



Catalog Number: 161142

---

Immunohistochemistry - ATF6 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using ATF6 antibody (40x lens).