

COG2 Polyclonal Antibody



Catalog Number: 161175

Product name

COG2 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to COG2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human COG2 (NP_031383.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:20 - 1:100

Validations



Immunofluorescence - COG2 Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using COG2 antibody . Blue: DAPI for nuclear staining.