

# FFAR3 Polyclonal Antibody



Catalog Number: 161181

## Product name

FFAR3 Polyclonal Antibody

## Specificity

Rat

## Antibody description

Polyclonal antibody to FFAR3

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 257-346 of human FFAR3 (NP\_005295.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



## Western blot - FFAR3 Polyclonal Antibody

Western blot analysis of extracts of rat skeletal muscle, using FFAR3 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 60s.