

# ACADL Polyclonal Antibody



Catalog Number: 161205

## Product name

ACADL Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to ACADL

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 31-210 of human ACADL (NP\_001599.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

## Validations




Western blot - ACADL Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ACADL antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.



Immunohistochemistry - ACADL Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver using ACADL antibody at dilution of



1:100 (40x lens).

# ACADL Polyclonal Antibody



Catalog Number: 161205

---

Immunofluorescence - ACADL Polyclonal Antibody

ACADL antibody at dilution of 1:100 (40x lens).

Blue: DAPI for nuclear staining.

Immunofluorescence analysis of HeLa cells using