

VEGFB Polyclonal Antibody



Catalog Number: 161213

Product name

VEGFB Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to VEGFB

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 25-207 of human VEGFB (NP_003368.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC, IF

Dilutions

WB 1:1000 - 1:2000
IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations



Western blot - VEGFB Polyclonal Antibody

Western blot analysis of extracts of U-87MG cells, using VEGFB antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



time: 90s.

Immunohistochemistry - VEGFB Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human mammary cancer using VEGFB antibody at dilution of 1:200 (40x lens).