

PHOSPHO2 Polyclonal Antibody



Catalog Number: 161279

Product name

Polyclonal

PHOSPHO2 Polyclonal Antibody

Ig Type**Specificity**

Rabbit IgG

Human

Applications**Antibody description**

WB

Polyclonal antibody to PHOSPHO2

Dilutions**Preparation**

WB 1:500 - 1:1000

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human PHOSPHO2 (NP_001008489.1).

Validations**Formulation**

Western blot - PHOSPHO2 Polyclonal Antibody

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Western blot analysis of extracts of various cell lines, using PHOSPHO2 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality