

TMEM126A Polyclonal Antibody



Catalog Number: 161286

Product name

TMEM126A Polyclonal Antibody

Specificity

Human, Rat

Antibody description

Polyclonal antibody to TMEM126A

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-195 of human TMEM126A (NP_115649.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - TMEM126A Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TMEM126A antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.