

# TLR8 Polyclonal Antibody



Catalog Number: 161348

## Product name

TLR8 Polyclonal Antibody

IHC 1:50 - 1:200

IF 1:50 - 1:200

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to TLR8

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 849-1041 of human TLR8 (NP\_619542.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB, IHC, IF

## Dilutions


WB 1:500 - 1:2000

## Validations



Western blot - TLR8 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TLR8 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.



Immunofluorescence - TLR8 Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using TLR8 antibody. Blue: DAPI for nuclear staining.