

AP3D1 Polyclonal Antibody



Catalog Number: 161422

Product name

AP3D1 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to AP3D1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 520-700 of human AP3D1 (NP_003929.4).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - AP3D1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using AP3D1 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.