

LRAT Polyclonal Antibody



Catalog Number: 161547

Product name

LRAT Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to LRAT

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-194 of human LRAT (NP_004735.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB 1:1000 - 1:2000

IF 1:50 - 1:200

Validations



Western blot - LRAT Polyclonal Antibody
Western blot analysis of extracts of various cell lines, using Lrat antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit .Exposure time: 30s.

Immunofluorescence - LRAT Polyclonal Antibody

Immunofluorescence analysis of MCF7 cells using LRAT antibody . Blue: DAPI for nuclear staining.