

ETV1 Polyclonal Antibody



Catalog Number: 161567

Product name

ETV1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to ETV1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 80-300 of human ETV1 (NP_004947.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC, IF

Dilutions

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations



Immunohistochemistry - ETV1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat pancreas using ETV1 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - ETV1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human breast cancer using ETV1 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - ETV1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse pancreas using ETV1 antibody at dilution

ETV1 Polyclonal Antibody



Catalog Number: 161567

of 1:100 (40x lens).



Immunofluorescence - ETV1 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using ETV1 antibody .