

MBTPS1 Polyclonal Antibody



Catalog Number: 161595

Product name

MBTPS1 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to MBTPS1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 803-1052 of human MBTPS1 (NP_003782.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC, IF

Dilutions

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations




Immunohistochemistry - MBTPS1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon cancer using MBTPS1 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - MBTPS1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon using MBTPS1 antibody at dilution of



1:100 (40x lens).


Immunohistochemistry - MBTPS1 Polyclonal Antibody

MBTPS1 Polyclonal Antibody



Catalog Number: 161595

Immunohistochemistry of paraffin-embedded human stomach using MBTPS1 antibody at



dilution of 1:100 (40x lens).

Immunofluorescence - MBTPS1 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using MBTPS1 antibody . Blue: DAPI for nuclear staining.