

# DNAJA3 Polyclonal Antibody



Catalog Number: 161600

## Product name

DNAJA3 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to DNAJA3

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 221-480 of human DNAJA3 (NP\_005138.3).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions


IF 1:50 - 1:200

## Validations



Immunofluorescence - DNAJA3 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using DNAJA3 antibody . Blue: DAPI for nuclear



staining.

Immunofluorescence - DNAJA3 Polyclonal Antibody

Immunofluorescence analysis of C6 cells using DNAJA3 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.