

CFHR3 Polyclonal Antibody



Catalog Number: 161609

Product name

CFHR3 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to CFHR3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 80-200 of human CFHR3 (NP_066303.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

Validations



Immunofluorescence - CFHR3 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using CFHR3 antibody . Blue: DAPI for nuclear staining.