

# UBIAD1 Polyclonal Antibody



Catalog Number: 161620

## Product name

UBIAD1 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to UBIAD1

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 200 to the C-terminus of human UBIAD1 (NP\_037451.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

## Validations



Western blot - UBIAD1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using UBIAD1 antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.



Immunohistochemistry - UBIAD1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using UBIAD1 antibody at dilution of 1:100 (40x lens).