

PHF11 Polyclonal Antibody



Catalog Number: 161626

Product name

PHF11 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to PHF11

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-292 of human PHF11 (NP_001035533.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations



Western blot - PHF11 Polyclonal Antibody

Western blot analysis of extracts of SW480 cells, using PHF11 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.




Immunohistochemistry - PHF11 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using PHF11 antibody at dilution of 1:100 (40x lens).

PHF11 Polyclonal Antibody



Catalog Number: 161626



Immunofluorescence - PHF11 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using PHF11 antibody .