

NAT8B Polyclonal Antibody



Catalog Number: 161628

Product name

NAT8B Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to NAT8B

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 75-165 of human NAT8B (NP_057431.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - NAT8B Polyclonal Antibody

Immunofluorescence analysis of A549 cells using NAT8B antibody . Blue: DAPI for nuclear staining.