

TET3 Polyclonal Antibody



Catalog Number: 161637

Product name

TET3 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to TET3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1556-1795 of human TET3 (NP_001274420.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC, IF

Dilutions

IHC 1:50 - 1:200

IF 1:50 - 1:200

Validations



Immunofluorescence - TET3 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using TET3 antibody . Blue: DAPI for nuclear staining.