

# GJA5 Polyclonal Antibody



Catalog Number: 161655

---

## Product name

GJA5 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to GJA5

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 229-358 of human GJA5 (NP\_859054.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions

IF 1:50 - 1:200

## Validations



Immunofluorescence - GJA5 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using GJA5 antibody . Blue: DAPI for nuclear staining.