

# PPP1CB Polyclonal Antibody



Catalog Number: 161673

## Product name

PPP1CB Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to PPP1CB

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-327 of human PPP1CB (NP\_002700.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

## Validations




Western blot - PPP1CB Polyclonal Antibody

Western blot analysis of extracts of SW480 cells, using PPP1CB antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry - PPP1CB Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human esophageal cancer using PPP1CB antibody



at dilution of 1:100 (40x lens).

# PPP1CB Polyclonal Antibody



Catalog Number: 161673

---

Immunofluorescence - PPP1CB Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using PPP1CB antibody . Blue: DAPI for nuclear staining.