

UPF1 Polyclonal Antibody



Catalog Number: 161680

Product name

UPF1 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to UPF1

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 270-370 of human UPF1 (NP_002902.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Immunohistochemistry - UPF1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver cancer using UPF1 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - UPF1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using UPF1 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - UPF1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse heart using UPF1 antibody at dilution of

UPF1 Polyclonal Antibody



Catalog Number: 161680

1:100 (40x lens).