

# ST3GAL4 Polyclonal Antibody



Catalog Number: 161682

## Product name

ST3GAL4 Polyclonal Antibody

## Specificity

Mouse, Rat

## Antibody description

Polyclonal antibody to ST3GAL4

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 27-267 of human ST3GAL4 (NP\_006269.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

## Validations



## Western blot - ST3GAL4 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ST3GAL4 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 60s.



## Immunohistochemistry - ST3GAL4 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat kidney using ST3GAL4 antibody at dilution of 1:200 (40x lens).