

NUTF2 Polyclonal Antibody



Catalog Number: 161696

Product name

NUTF2 Polyclonal Antibody

Specificity

Human, Rat

Antibody description

Polyclonal antibody to NUTF2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-127 of human NUTF2 (NP_005787.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

Validations



Immunofluorescence - NUTF2 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using NUTF2 antibody . Blue: DAPI for nuclear staining.