

# TIMM17A Polyclonal Antibody



Catalog Number: 161697

## Product name

TIMM17A Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to TIMM17A

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-171 of human TIMM17A (NP\_006326.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

## Validations




Western blot - TIMM17A Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TIMM17A antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 30s.



Immunohistochemistry - TIMM17A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver cancer using TIMM17A antibody at



dilution of 1:100 (40x lens).

# TIMM17A Polyclonal Antibody



Catalog Number: 161697

---

Immunofluorescence - TIMM17A Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using TIMM17A antibody . Blue: DAPI for nuclear staining.