

Eif2ak2 Polyclonal Antibody



Catalog Number: 161700

Product name

Eif2ak2 Polyclonal Antibody

Specificity

Human, Rat

Antibody description

Polyclonal antibody to Eif2ak2

Preparation

Antigen: Recombinant protein of mouse Eif2ak2

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:200

Validations



Immunohistochemistry - Eif2ak2 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat brain using Eif2ak2 antibody at dilution of 1:200



(40x lens).

Immunohistochemistry - Eif2ak2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung cancer using Eif2ak2 antibody at dilution of 1:200 (40x lens).



Immunohistochemistry - Eif2ak2 Polyclonal Antibody

Eif2ak2 Polyclonal Antibody



Catalog Number: 161700

Immunohistochemistry of paraffin-embedded

human liver cancer using Eif2ak2 antibody at dilution of 1:200 (40x lens).