

# P4HTM Polyclonal Antibody



Catalog Number: 161705

---

## Product name

P4HTM Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to P4HTM

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 294-563 of human P4HTM (NP\_808807.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions

IF 1:50 - 1:200

## Validations



Immunofluorescence - P4HTM Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using P4HTM antibody . Blue: DAPI for nuclear staining.