

SPDYA Polyclonal Antibody



Catalog Number: 161714

Product name

SPDYA Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to SPDYA

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-313 of human SPDYA (NP_877433.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

Validations



Western blot - SPDYA Polyclonal Antibody

Western blot analysis of extracts of A549 cells, using SPDYA antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 10s.



Immunofluorescence - SPDYA Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using SPDYA antibody. Blue: DAPI for nuclear staining.