

ZIC3 Polyclonal Antibody



Catalog Number: 161737

Product name

ZIC3 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to ZIC3

Preparation

Antigen: A synthetic peptide of human ZIC3

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, FC


Dilutions

WB 1:500 - 1:1000

IHC 1:50 - 1:100

FC 1:20 - 1:50

Validations



Western blot - ZIC3 Polyclonal Antibody

Western blot analysis of extracts of HL-60 cells, using ZIC3 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.