

DEFB132 Polyclonal Antibody



Catalog Number: 161744

Product name

DEFB132 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to DEFB132

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human DEFB132 (NP_997352.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Immunohistochemistry - DEFB132 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human colon using DEFB132 antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - DEFB132 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human esophagus using DEFB132 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - DEFB132 Polyclonal Antibody

DEFB132 Polyclonal Antibody



Catalog Number: 161744

Immunohistochemistry of paraffin-embedded

mouse kidney using DEFB132 antibody at dilution of 1:100 (40x lens).