

MAOA Polyclonal Antibody



Catalog Number: 161893

Product name

MAOA Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to MAOA

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human MAOA (NP_000231.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

Validations



Immunofluorescence - MAOA Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using MAOA antibody . Blue: DAPI for nuclear staining.