

# MSH4 Polyclonal Antibody



Catalog Number: 161895

## Product name

MSH4 Polyclonal Antibody

## Specificity

Human, Rat

## Antibody description

Polyclonal antibody to MSH4

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 697-936 of human MSH4 (NP\_002431.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

## Validations



Western blot - MSH4 Polyclonal Antibody

Western blot analysis of extracts of SKOV3 cells, using MSH4 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.



Immunofluorescence - MSH4 Polyclonal Antibody

Immunofluorescence analysis of rat testis using MSH4 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.